



Liquid Handling

ADDRESS

NO.9, SIDDHIVINAYAK COMPLEX, KHANIVALI, TAL: WADA, DIST: PALGHAR, MAHARASHTRA INDIA - 421312



WWW.SMARTLABINN.COM


SMARTLABINNOVATORS@GMAIL.COM



+91-97 6962 4321

DLAB Liquid Handling range includes Micropipettes (Electronic Pipettes & Mechanical Pipettes), Pipette Filler, Pipette Controller and Stepper.

List of DLAB models:

Series	Image	Model	Application	Operation Method	Volume Range	Channel Option	Autoclavable	Page
Mechanical Pipette		HiPette	Aqueous solution	Manual	0.1-10000µL	Single channel	YES	p.5
		MicroPette	Aqueous solution	Manual	0.1-10000µL	Single Channel Adj. & Fixed;	YES (lower part)	p.8
					0.5-300µL	8 Channel Adj.; 12 Channel Adj.;		
		MicroPette Plus	Aqueous solution	Manual	0.1-5000µL	Single Channel Adj. & Fixed;	YES	p.7
					0.5-300µL	8 Channel Adj.; 12 Channel Adj.;		
		MiniPette	Aqueous solution	Manual	25-200µL	Single Channel Fixed	NO	p.18
Electronic Pipette		dPette+	Aqueous solution	Electronic	0.5-1000µL	Single Channel Adj.	YES (lower part)	p.11
					0.5-300µL	8 Channel Adj.		
		dPette	Aqueous solution	Electronic	0.5-1000µL	Single Channel Adj.	YES (lower part)	p.13
Pipette Filler		Levo Plus	Aqueous solution	Electronic	0.1-100mL	Single channel	YES (silicon adapter&nozzle)	p.15
		Levo Me	Aqueous solution	Electronic	0.1-100mL	Single channel	YES (silicon adapter&nozzle)	p.17
Pipette Controller		Levo	Aqueous solution	Manual	0.1-100mL	Single channel	YES (silicon adapter&nozzle)	p.18
Pipette Pump		Levo E	Aqueous solution	Manual	2-25mL	Single channel	NO	p.18
Stepper		StepMate	Aqueous solution, high viscosity or volatile solution	Manual	0.5-50mL	Single channel	NO	p.19

In addition to Pipettes, Liquid Handling series also include various Dispensers, Electronic Burette and Vacuum Aspiration System, which are used for the pipetting, dispensing, and collection of reagents with a wider range.

List of DLAB models:



Product Name	Bottle-Top Dispenser	New Bottle-Top Dispenser	Electronic Bottle-Top Dispenser	Electronic Burette	Vacuum Aspiration System	Economical Vacuum Aspirator
Model	DispensMate	DispensMate-Pro	dFlow	dTrite	SafeVac	EcoVac
Application	Reagent pipetting and liquid dispensing			Used in titration operations in chemical analysis, food industry, environmental analysis	Liquid collection and storage treatment, including cell culture medium, suspension, supernatant.	
Operation	Manual	Manual	Electronic	Electronic	Electronic	Electronic
Volume Range	0,5-5mL 1,0-10mL 2,5-25mL 5,0-50mL	0,5-5mL 1,0-10mL 2,5-25mL 5,0-50mL 10-100mL	0,1-99,9mL	0,01-99,99mL	Vacuum range: 0-600mbar Bottle volume:4L	Vacuum range: 0-500mbar Bottle volume:2L
Reagent Recovery Function	NO	Yes	Yes	Yes	NO	NO
Features	-	-Pro Glass piston (Part volume)	-	-	Sensitive level sensor, Self-closing connector	1L bottle optional, Switch manual and continuous aspiration
Fully autoclavable	The part where the liquid flows	Whole	-	-	The part where the liquid flows	The part where the liquid flows
Accessories	Screw adapter	Screw adapter	Spiral dispensing pipe	Spiral dispensing pipe	Handle adapter	Handle adapter
Page	P.24	P.23	P.22	P.21	P.25	P.27

HiPette

Fully Autoclavable Mechanical Pipette

The DLAB new generation fully autoclavable mechanical pipette has the best ergonomic design and ultralight weight. It has been designed keeping in mind the human form and human factors to offer smooth and effortless ergonomic pipetting experience. It is UV sterilisable and resistant to strong chemical corrosion. In addition, the volume lock helps protect from volume changing accidentally, and achieve reliable pipetting of smallest volume of liquids.



Ultra-light ergonomic design

Light weight and light pipetting force with soft range adjustment and piston movement offers effortless & fatigue-less pipetting experience.



Accurate pipetting

Special material and the updated components & technologies effectively reduces the operating force for accuracy, precision and reliability.



Complete sterilization

The best in class pipette material used is fully autoclavable. It gives better protection against chemical & physical corrosion, and offers enhanced UV resistance for higher durability.



Specifications

SINGLE-CHANNEL ADJUSTABLE VOLUME PIPETTES
(9 MODELS, COVERING 0.1 μ L~10mL.)

Volume Range	Increment	Test Volume	Accuracy Error		Precision Error	
μ L	μ L	μ L	μ L	%	μ L	%
0.1~2.5	0.002	2.5	± 0.04	± 1.4	± 0.02	± 0.7
		1.25	± 0.03	± 2.5	± 0.02	± 1.5
		0.25	± 0.03	± 12.0	± 0.02	± 6.0
0.5~10	0.01	10	± 0.10	± 1.0	± 0.04	± 0.4
		5	± 0.08	± 1.5	± 0.04	± 0.8
		1	± 0.03	± 2.5	± 0.02	± 1.8
2~20	0.02	20	± 0.20	± 1.0	± 0.06	± 0.3
		10	± 0.12	± 1.2	± 0.06	± 0.6
		2	± 0.10	± 5.0	± 0.03	± 1.5
10~100	0.1	100	± 0.80	± 0.8	± 0.20	± 0.2
		50	± 0.50	± 1.0	± 0.15	± 0.3
		10	± 0.30	± 3.0	± 0.10	± 1.0
20~200	0.2	200	± 1.20	± 0.6	± 0.40	± 0.2
		100	± 1.00	± 1.0	± 0.30	± 0.3
		20	± 0.50	± 2.5	± 0.14	± 0.7
30~300	0.2	300	± 1.80	± 0.6	± 0.60	± 0.2
		150	± 1.50	± 1.0	± 0.45	± 0.3
		30	± 0.75	± 2.5	± 0.21	± 0.7
100~1000	1	1000	± 6.00	± 0.6	± 2.00	± 0.2
		500	± 5.00	± 1	± 1.00	± 0.2
		100	± 3.00	± 3	± 0.60	± 0.6
500~5000**	5	5000	± 30.00	± 0.6	± 10.00	± 0.2
		2500	± 15.00	± 0.6	± 7.50	± 0.3
		500	± 12.00	± 1.2	± 3.00	± 0.3
1000~10000**	10	10000	± 60.00	± 0.6	± 20.00	± 0.2
		5000	± 40.00	± 0.8	± 10.00	± 0.2
		1000	± 30.00	± 1.5	± 6.00	± 0.3

** Coming soon

User calibration should refer to the industrial standard ISO8655-2.


MicroPette Plus

121°C
FULLY
Autoclavable



MicroPette

121°C
Half
Autoclavable



Features

- Available as Fixed Volume and Variable Volume formats
- Streamlined shape design
- Comprehensive volume range from 0.1µL to 10mL
- Easy calibration and maintenance
- 8 and 12 channel pipette options
- Dispensing head rotates for effortless pipetting convenience
- Individual piston and tip cone assembly
- Spring loaded tip cones for easy cleaning and maintenance
- Compound material-made tip cone secures high sealing performance



Mechanical Pipette Volume Selection

Specifications

This volume list is for MicroPette and MicroPette plus (Adjustable and Fixed volume)

Single-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume	Error limits in accordance with ISO8655-2			
			Accuracy error		Precision error	
μL	μL	μL	μL	%	μL	%
0.1–2.5	0.05	2.5	± 0.0625	± 2.50	± 0.05	± 2.00
		1.25	± 0.0375	± 3.00	± 0.0375	± 3.00
		0.25	± 0.03	± 12.00	± 0.015	± 6.00
0.5–10	0.1	10	± 0.1	± 1.00	± 0.08	± 0.80
		5	± 0.075	± 1.50	± 0.075	± 1.50
		1	± 0.025	± 2.50	± 0.015	± 1.50
2–20	0.5	20	± 0.18	± 0.09	± 0.08	± 0.04
		10	± 0.12	± 1.20	± 0.1	± 1.00
		2	± 0.06	± 3.00	± 0.04	± 2.00
5–50	0.5	50	± 0.3	± 0.60	± 0.15	± 0.30
		25	± 0.225	± 0.90	± 0.15	± 0.60
		5	± 0.1	± 2.00	± 0.1	± 2.00
10–100	1	100	± 0.8	± 0.80	± 0.15	± 0.15
		50	± 0.5	± 1.00	± 0.2	± 0.40
		10	± 0.3	± 3.00	± 0.15	± 1.50
20–200	1	200	± 1.2	± 0.60	± 0.3	± 0.15
		100	± 0.8	± 0.80	± 0.3	± 0.30
		20	± 0.6	± 3.00	± 0.2	± 1.00
50–200	1	200	± 1.2	± 0.60	± 0.3	± 0.15
		100	± 0.8	± 0.80	± 0.3	± 0.30
		50	± 0.5	± 1.00	± 0.2	± 0.40
100–1000	5	1000	± 6	± 0.60	± 2	± 0.20
		500	± 3.5	± 0.70	± 1.25	± 0.25
		100	± 2	± 2.00	± 0.7	± 0.70
200–1000	5	1000	± 6	± 0.60	± 2	± 0.20
		500	± 3.5	± 0.70	± 1.25	± 0.25
		200	± 1.8	± 0.90	± 0.6	± 0.30
1000–5000	50	5000	± 25	± 0.50	± 7.5	± 0.15
		2500	± 15	± 0.60	± 7.5	± 0.30
		1000	± 7	± 0.70	± 3	± 0.30
2000–10000	100	10000	± 60	± 0.60	± 20	± 0.20
		5000	± 60	± 1.20	± 15	± 0.30
		2000	± 60	± 3.00	± 12	± 0.60

User calibration should refer to the industrial standard ISO8655-2.

Fixed Volume Pipettes						
Volume Range	Increment	Test Volume	Accuracy error		Precision error	
μL	–	μL	μL	%	μL	%
5	–	5	± 0.065	± 1.3	± 0.06	± 1.2
10	–	10	± 0.08	± 0.8	± 0.08	± 0.8
20	–	20	± 0.12	± 0.6	± 0.1	± 0.5
25	–	25	± 0.125	± 0.5	± 0.075	± 0.3
50	–	50	± 0.25	± 0.5	± 0.15	± 0.3
100	–	100	± 0.5	± 0.5	± 0.3	± 0.3
200	–	200	± 0.8	± 0.4	± 0.4	± 0.2
250	–	250	± 1.0	± 0.4	± 0.5	± 0.2
500	–	500	± 1.5	± 0.3	± 1.0	± 0.2
1000	–	1000	± 3.0	± 0.3	± 2.0	± 0.2
2000	–	2000	± 6.0	± 0.3	± 3.0	± 0.15
5000	–	5000	± 15	± 0.3	± 7.5	± 0.15

8-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume	Accuracy error		Precision error	
μL	μL	μL	μL	%	μL	%
0.5–10	0.1	10	± 0.15	± 1.50	± 0.15	± 1.50
		5	± 0.125	± 2.50	± 0.125	± 2.50
		1	± 0.04	± 4.00	± 0.04	± 4.00
5–50	0.5	50	± 0.5	± 1.00	± 0.25	± 0.50
		25	± 0.375	± 1.50	± 0.25	± 1.00
		5	± 0.15	± 3.00	± 0.1	± 2.00
50–300	5	300	± 2.1	± 0.70	± 0.75	± 0.25
		150	± 1.5	± 1.00	± 0.75	± 0.50
		50	± 0.75	± 1.50	± 0.4	± 0.80

12-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume	Accuracy error		Precision error	
μL	μL	μL	μL	%	μL	%
0.5–10	0.1	10	± 0.15	± 1.50	± 0.15	± 1.50
		5	± 0.125	± 2.50	± 0.125	± 2.50
		1	± 0.04	± 4.00	± 0.04	± 4.00
5–50	0.5	50	± 0.5	± 1.00	± 0.25	± 0.50
		25	± 0.375	± 1.50	± 0.25	± 1.00
		5	± 0.15	± 3.00	± 0.1	± 2.00
50–300	5	300	± 2.1	± 0.70	± 0.75	± 0.25
		150	± 1.5	± 1.00	± 0.75	± 0.50
		50	± 0.75	± 1.50	± 0.4	± 0.80

* DLAB specifications are used as guidelines and the user calibration should refer to the industrial standard ISO 8655

dPette+

Multi-functional 8-channel Electronic Pipette

DLAB multi-functional high performance 8-channel electronic pipette offers productive pipetting with its easy to understand operations. Its uniquely lightweight streamlined design ensures an effortless transfer of multiple samples with an increased throughput and data reproducibility.



• Easy Operation

Intuitive interface for setting functions and parameters.

• Ergonomic Design

Low operation forces for complete work, which is exceptionally fatigue-free.

• Convenient and Versatile

360° pipetting. Double knobs for simple and versatile control. Easy loading tip cones offers smooth and leak free operation.

• Dual charging modes

Use of USB charger or the charging stand to ensure uninterrupted use.

Three easy steps to operate dPette+

1

Long press the Parameter Knob for 2 seconds to Start

2


Quickly turn the Function Wheel to Activate, Rotate to Switch Pipetting, Continuous Dispensing, and other Function Settings

3

After quickly turning the Parameter Knob to Unlock, Turn to Adjust the Volume, Press to Pipette, Long press to Mix


Pipetting function

Parameter Knob press → Pipetting,
Long press → Mixing




Dispensing function

Dispensing times = Pipette max volume /
Single sample volume



Other function

Pipetting speed adjustment
Key tone adjustment



Complimentary Pipette Holder

Supports dual charging mode



USB Charging

Automatic Calibration Support
Connect it to the computer,
with the free software
provided for calibration



Contact Charging



Specifications

Channels	Volume Range	Increment	Test Volume	Accuracy error		Precision error	
	μL	μL	μL	μL	%	μL	%
8	0.5-10	0.01	10	± 2.00	± 0.20	± 1.0	± 0.1
			1	± 8.0	± 0.08	± 5.00	± 0.05
8	10-100	0.1	100	± 0.80	± 0.80	± 0.30	± 0.3
			10	± 3.00	± 0.3	± 2.00	± 0.2
8	30-300	1	300	± 0.60	± 1.80	± 0.3	± 0.9
			30	± 3.00	± 0.90	± 1.0	± 0.3

* SD=Standard Deviation * CV=Co-efficient of Variation

* DLAB specifications are used as guidelines and the user calibration should refer to the industrial standard ISO 8655.

dPette is an innovative Electronic Pipette developed by DLAB. It combines manual pipette features, as ergonomics and lightweight, with electronic pipette features, such as labor-saving and high accuracy, offering new pipetting experience to users.

Parameter Knob

Function Wheel

Display

ELECTRONIC PIPETTE

Pipette Tip cone

Pipetting Mode
Dispensing Mode
Diluting Mode

Buzzer & Settings

Battery Indicator

Volume Range
Parameter

Unit

Indicator

Power off



dPette+

Multi functional Electronic Pipette

dPette

Simple Electronic Pipette

Features

- Motor driven digitally control pipette with multifunctions
- Easy Operations and Handling
- Intuitive menu interface settings; functions and operations
- Minimal force for pipetting operation
- High performance ensuring accuracy & repeatability
- 2 buttons for all operational settings
- Adjustable speed for aspiration and dispensing
- Li-ion battery enable longer operation time
- Convenient and Versatile
- Self-calibration applicable



dPette+
For Pipetting,
Mixing, Stepper
and Dilution



dPette
For Pipetting
and Mixing

Specifications

Channels	Volume Range	Increment	Test Volume	accuracy error		precision error	
	µL	µL	µL	µL	%	µL	%
1	0.5–10	0.01	10	±0.10	±1.00	±0.05	±0.50
			1	±0.035	±3.50	±0.03	±3.00
1	5–50	0.1	50	±0.40	±0.80	±0.15	±0.30
			5	±0.15	±3.00	±0.125	±2.50
1	30–300	1	300	±1.80	±0.60	±0.60	±0.20
			30	±0.09	±3.00	±0.21	±0.70
1	100–1000	5	1000	±6.00	±0.60	±2.00	±0.20
			100	±3.00	±3.00	±0.60	±0.60

DLAB specifications are used as guidelines and the user calibration should refer to the industrial standard ISO 8655

Levo Plus

Pipette Filler

Ergonomic fast pipetting solution with accurate speed control



Features

- Intuitive and convenient speed with 6 levels
- Lightweight pipette filler compatible with glass or plastic serological pipettes (0.1 – 100 mL)
- Ergonomic design for fatigue-free and effective streamlined pipetting
- LCD screen displays battery status and pipetting speed
- Quiet and powerful motor and pump fills 25 mL pipette in less than 5 seconds
- Long-lasting lithium-ion battery



Select pipetting speed by pressing +/- button

Large LCD display for low battery warning and speed settings

Aspiration button

Operating button for motor and gravity dispensing

Hydrophobic filter 0.45µm

Silicone adapter for 1.0–100mL pipette tip, fully autoclavable



Specifications

Aspiration Speed	25mL<5s (6 Level)
Dispensing Speed	Motor (6 shift)/Gravity
Battery	Lithium-ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2–3 Hours
Pipette Types	Glass or Plastic Pipette (0.1–100mL), Pasteur Pipettes
Filter	0.45µm Hydrophobic
Weight	200g

Levo ME

Pipette Filler

Elegant design & reliable performance,
with long-life battery

Features

- Enables single-handed operation with minimum effort
- Speed adjustment simply with finger tips
- Individually set aspiration (upper button) and dispense (lower button) speeds
- Compatible with glass or plastic serological pipettes (0.1–100mL)
- 0.45µm replaceable hydrophobic filter
- High capacity Li-ion battery enable long operation time

Specifications

Aspiration Speed	25ml<7s
Dispensing Speed	Motor/Gravity
Battery	Lithium-ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2–3 Hours
Pipette Types	Glass or Plastic Pipette (0.1–100mL), Pasteur Pipettes
Filter	0.45µm Hydrophobic
Weight	200g




Levo

Pipette Controller

Features

- Reduce hands-on time
- Compatible with glass or plastic serological pipettes (0.1 - 100mL)
- 3µm replaceable hydrophobic filter
- Easy to clean and maintain



Mini Pipette

Features

- Ergonomic design provides excellent operation experience;
- Smart body, ideal to combine with IVD reagents;
- Compatible with universal 200µL tips;

Volumes
200µL 100µL 75µL 50µL 30µL 25µL OEM

Levo E

Pipette Pump

Designed for single-handed operation, Pipette Pump identifies itself with a thumbwheel, allowing precise aspiration or dispensing. The feature of side lever enables complete liquid dispensing, without a droplet remaining.



Features

- Volume capacity 2mL, 10mL and 25mL
- Colour-coded by volume with green, blue and red
- Thumb wheel guarantees precise operation
- Resistant to acids and alkalis
- Easy for cleaning and maintenance

StepMate

Stepper

Positive-displacement, continuously adjustable, repetitive dispenser



Features

- Slim, light and robust
- wide dispensing from 10µL to 5mL
- Up to 48 dispensing steps in succession
- Equipped with durable tip insertion lever
- With dispensing tips volume sizes from 0.5mL to 50mL

Applications

Ideal for any application requiring repetitive dispensing in the clinical, hospital, biological, chemical, food and beverage, pharmaceutical industry.

How to select dispensing volume

Find the required dispensing volume in the table.

Use the volume selection dial to set the required volume.

Select, and insert a dispensing tip.

Follow the maximum dispensing steps as stated.

DLAB dispensing tips are convenient for viscous or foaming solutions, volatile liquids and applications where repetitive dispensing of volumes in long series is required

Dispensing tips are made from high quality materials without soft lubricants or other additives.

Specifications

Note: The first and last dispensing liquid should be discarded for precision.



Dial setting		1	2	3	4	5			
Number of steps		48	23	15	11	8			
Cat. No.	Dispensing tips volume	Sample volume					Accuracy error	Precision error	Description
	μL	μL					%	%	
17900032	0.50	10	20	30	40	50	± 0.8	≤ 0.7	0.5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900033	1.25	25	50	75	100	125	± 0.8	≤ 0.5	1.25mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900034	2.50	50	100	150	200	250	± 0.8	≤ 0.4	2.5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900035	5.00	100	200	300	400	500	± 0.4	≤ 0.3	5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900036	12.50	250	500	750	1000	1250	± 0.3	≤ 0.3	12.5mL, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900037	25.00	500	1000	1500	2000	2500	± 0.2	≤ 0.3	25mL, piston PE-HD, cylinder PP, 1 packing, 2pcs/packing with 1 pc adapter
17900038	50.00	1000	2000	3000	4000	5000	± 0.2	≤ 0.2	50mL, piston PE-HD, cylinder PP, 1 packing, 2pcs/packing with 1 pc adapter