

Online Monitoring Systems



Made in India



Contact Us

smartlabinn@gmail.com

info@smartlabinn.com



**ISO 9001:2015
ISO 13485:2016
Certified Organization**

**We are one of the eminent manufacturers, supplier, service providers of
Online Monitoring Systems,
Medical, Pharmaceutical and Laboratory Instruments & Equipments.**

**We always try to keep us in touch with the updated technologies in the
field of Scientific Instruments to provide the best technology for our
valuable customers.**

**We continue to grow, adding new products based on a strong foundation
of values, technical innovation, high quality, customer satisfaction and of
course our perseverance.**


**Contact Us
smartlabinn@gmail.com
info@smartlabinn.com**

Online Temperature Monitoring System



CE complaint
RoHS complaint

Temperature Range	-200 to +700 Deg C
High/Low Alarms	Audio & Visual
Sensor Type	RTD PT 100
On unit Display	Yes
Web Dashboard	Yes
Data Logs	Yes - Online
Device Logs	Yes - Online
Connectivity	Wifi / Internal Server
Login Access	1 Admin & 5 users
Alerts	email, sms, whatsapp
Battery Backup	Yes - Optional




Contact Us
smartlabinn@gmail.com
info@smartlabinn.com

Online Humidity-Temperature Monitoring System



Humidity Range	0 to 100 % RH
Temperature Range	-40 to +125 Deg C
High/Low Alarms	Audio & Visual
On unit Display	Yes
Web Dashboard	Yes
Data Logs	Yes - Online
Device Logs	Yes - Online
Connectivity	Wifi / Internal Server
Login Access	1 Admin & 5 users
Alerts	email, sms, whatsapp
Battery Backup	Yes - Optional




Contact Us
smartlabinn@gmail.com
info@smartlabinn.com

Online Room Oxygen Monitoring System



CE complaint
RoHS complaint

Range	0 to 25 % vol
High/Low Alarms	Audio & Visual
On unit Display	Yes
Web Dashboard	Yes - Optional
Data Logs	Yes - Optional
Device Logs	Yes - Optional
Connectivity	Wifi / Internal Server
Login Access	Optional
Alerts	email,sms - Optional
Battery Backup	Yes - Optional



Contact Us
smartlabinn@gmail.com
info@smartlabinn.com




Available Modules

- Temperature (Range -200 to +700 Deg celsius)
- Relative Humidity (Range 0 to 100%)
- Oxygen Level (Range 0% to 25% & upto 100%)
- CO2 Level (Range 0 to 0.5%, 0 to 20%)
- Line Pressure Monitors (Available in bar, psi, etc)
- Liquid Nitrogen Low Level Alarm


We also provide customized services as per customer requirement.

Features




- Remote Monitoring
- Realtime Monitoring
- Dedicated Display on each module
- Wireless Technology
- Cloud as well as Internal server option
- High & Low set able alarms
- 21 CFR Part 11 Complaint
- Secured login access
- Data Logs
- Device Logs
- Graphical representation
- Email, SMS, Whatsapp alerts

Internal Server Based Model




SMARTLAB

INNOVATORS



Dash Board



Data Logs

SMARTLAB INNOVATORS www.smartlabinn.com -80 Degree 4 OMS10056TMP				
OMS10056TMP	24-03-2023	02:20 PM	66.0 °C	--
OMS10056TMP	24-03-2023	02:20 PM	66.0 °C	--
OMS10056TMP	24-03-2023	02:22 PM	67.9 °C	--
OMS10056TMP	24-03-2023	02:26 PM	67.8 °C	--
OMS10056TMP	24-03-2023	02:23 PM	67.9 °C	--
OMS10056TMP	24-03-2023	02:24 PM	68.3 °C	--
OMS10056TMP	24-03-2023	02:23 PM	68.7 °C	--
OMS10056TMP	24-03-2023	02:22 PM	70.0 °C	--
OMS10056TMP	24-03-2023	02:15 PM	-76.9 °C	--
OMS10056TMP	24-03-2023	02:00 PM	-77.6 °C	--
OMS10056TMP	24-03-2023	01:45 PM	-76.9 °C	--
OMS10056TMP	24-03-2023	01:30 PM	-78.5 °C	--
OMS10056TMP	24-03-2023	01:15 PM	-78.2 °C	--
OMS10056TMP	24-03-2023	01:00 PM	-77.9 °C	--
OMS10056TMP	24-03-2023	12:45 PM	-77.5 °C	--
OMS10056TMP	24-03-2023	12:30 PM	-76.9 °C	--
OMS10056TMP	24-03-2023	12:15 PM	-77.0 °C	--
OMS10056TMP	24-03-2023	11:45 AM	-79.7 °C	--
OMS10056TMP	24-03-2023	11:30 AM	-79.8 °C	--
OMS10056TMP	24-03-2023	11:15 AM	-79.8 °C	--
OMS10056TMP	24-03-2023	11:00 AM	-79.6 °C	--
OMS10056TMP	24-03-2023	10:45 AM	-79.8 °C	--
OMS10056TMP	24-03-2023	10:30 AM	-79.7 °C	--
OMS10056TMP	24-03-2023	10:15 AM	-79.7 °C	--
OMS10056TMP	24-03-2023	10:00 AM	-79.7 °C	--

Our Installations



Contact Us
smartlabinn@gmail.com
info@smartlabinn.com



Our Other Services

Laboratory & Pharma Instruments

- **AMC Services**
- **Service and Repairs**
- **Calibration / Validation**
- **Customized Instruments**
- **Preventive Maintenance Services**

Contact Us
smartlabinn@gmail.com
info@smartlabinn.com



Our Other Instruments

- Biosafety Cabinets
- BOD Incubators
- CO2 Incubators
- Water Baths
- Centrifuges
- Deep Freezers
- Ultra Low Temperature Freezers
- Vortex
- Rotators
- Mini centrifuges
- Hot Plates
- Dry Bath
- Micro Pipettes
- Probe Sonicators

Contact Us
smartlabinn@gmail.com
info@smartlabinn.com

Application



Laboratories



Biorepositories



Cryogenics



Blood Banks



IVF



Hospitals



Cold Rooms



Server Rooms



Industries



Warehouses



Mining

Contact Us
smartlabinn@gmail.com
info@smartlabinn.com



Our Esteemed Clients

Anuh Pharma AIIMS Nagpur NCCS

Pune University MITCON Serum Institute of India

Cryoviva Singapore NCL-Venture Center Pensieve HealthCare

and others.....



Address:
No.9, Siddhivinayak Complex,
Khaniwali, Tal: Wada, Dist: Palghar
Maharashtra, India 401204

+91-97 6962 4321

Contact Us
www.smartlabinn.com
smartlabinn@gmail.com
info@smartlabinn.com